

Meywalk Gait Trainer
Order Form 2008

Pacific Rehab, Inc
36805 N Nevermind Trail, Box 5406
Carefree AZ, 85377-5406

Toll Free: 888.222.9040
Fax: 480-575-7907

Account Information

Company Name: _____
Address: _____
City, State Zip: _____
Phone: _____
Fax: _____
Email: _____
Purchase Order: _____

Mark For: _____
Ship To: _____
Address: _____
State Zip: _____
Phone: _____

Order Information

To insure proper gait trainer sizing, please include required measurements below

Miniwalk I Gait Trainer
Seat Height: 10 1/2" - 19"
Trunk Support Height: 20" - 29"
Ring Size: 8 1/4" w x 8" d
Weight Capacity: 90 lbs

106061 Miniwalk I Gait Trainer 2,475.00 Ea
Includes Castor Swivel Locks
Color: Turquoise

Miniwalk II Gait Trainer
Seat Height: 14 - 25"
Trunk Support Height: 27-35"
Ring Size: 8 1/4" w x 8" d
Weight Capacity: 90 lbs

106062 Miniwalk II Gait Trainer 2,475.00 Ea
Includes Castor Swivel Locks
Color: Turquoise

Miniwalk I Accessories

106070 Antitip support 225.00 Pr
106080 Legguides 285.00 Pr
106075 Leg Separation Panel 450.00 Ea
106085 Hip Pads 475.00 Pr
106089 Non-reverse brakes 275.00 Pr
106096 Drag brakes 224.00 Pr

Meywalk 2000 Gait Trainer Medium
Seat Height: 21"-29"
Trunk Support Height: 33"-40"
Ring Size: 12 3/4" w x 11 1/2" d
Weight Capacity: 220lbs


104175 Meywalk 2000 Medium 2,997.00 Ea
Color: Blue

Meywalk 2000 Gait Trainer Large
Seat Height: 27" -37"
Trunk Support Height: 41" - 51"
Ring Size: 14 1/2" w x 13 3/4" d
Weight Capacity: 220lbs

104173 Meywalk 2000 Large 2,997.00 Ea
Color: Blue

Meywalk Accessories

104023 Non-reverse brakes 275.00 Pr
104044 Leg separation plate 450.00 Ea
104110 Height reducing fittings, -3" 115.00 Pr
104142 Extension Tubes, +2.3", Large Only 30.00 Pr
104166 Drag Brakes 224.00 Pr
104167 Antitip support 225.00 Pr
104168 Legguides 285.00 Pr
104215 Hip Pads 475.00 Pr
104997 Front Castor Swivel Locks 100.00 Pr



Required Measurements:

- A: Ring Height** _____
Measure user from floor to axilla to indicate ring height.
- B: Inseam Height:** _____
Measure user from floor to inseam, with user standing. Indicates seat height.
- C: Trunk Width:** _____
Measure user across trunk, at axilla height, side to side in a straight line, to indicate proper ring size.